



General catalogue Operating chairs and stools

Ergonomic seating comfort
and concentrated work

Welcome to BRUMABA

Reliable quality for your health

As a family-owned company, we specialize in the development and production of ergonomic surgical tables, chairs, and stools. By providing top-quality and innovative products, we aim to improve the health and well-being of medical professionals and patients worldwide.

At BRUMABA, we value sustainability: Our products have been manufactured in our Bavarian home in Southern Germany since 1980. We use only the highest quality materials and reliable technology, ensuring our products have a long service life.

Experience our leading edge in ergonomics, comfort and efficiency for yourself.

Reliability that counts.

Stability that lasts.

BRUMABA - your partner in quality.

Operating stool
COMFORT



MADE IN GERMANY



SINCE 1980



Operating chair
BALANCE



Operating chair
BALANCE ADVANCE

BRUMABA Operating Chairs and Stools: An Investment in Your Health

For surgeons and healthcare professionals, maintaining an ergonomic sitting posture is essential not only for comfort but also for professional efficiency. Especially during long and demanding procedures, proper ergonomics promote better blood circulation, enhanced concentration, and reduced muscle tension and fatigue.

BRUMABA surgical stools and chairs are designed to adapt flexibly to individual body measurements, serving as an ideal complement to our operating tables. Their ergonomic design promotes a back-friendly posture, minimizing physical strain during surgeries and assisting roles—a sustainable investment in both your personal health and that of your team.

Your benefits:

- Optimal back support with infinitely adjustable back sections
- Precisely adjustable armrests for optimal support during surgeries and treatments
- Patented ergonomic seat surfaces for consistently secure and comfortable seating
- Individually adjustable seat tilt angle for healthy blood circulation in the legs
- Maintaining a sterile working environment with foot-adjustable seat height

Investing in BRUMABA operating chairs and stools ensures not only comfort but also supports your long-term health and professional performance.





Optimal
back support

Multi-functional
armrest

Ergonomic
seat

Vein-friendly
sitting

Operating chair **BALANCE**

Achieving perfect balance: Form, function, and ergonomics

BALANCE F

Standout features

- Body-adaptive seat and backrest for ergonomic support
- High-quality cushions for comfortable and healthy sitting
- The ergonomic seat design promotes healthy blood circulation in your legs, enabling you to sit comfortably for extended periods of time
- The cushions are seamlessly designed and the wheels are enclosed, allowing effortless cleaning and ensuring impeccable hygiene
- Optional column elevation (100 mm) (ART. NO. E.001558)



Optional with brake



Height adjustment with foot release, ideal for sterile work

Technical data

Min. height	500 mm
Max. height	695 mm
Chair base diameter	610 mm
Seat cushion length	460 mm
Seat cushion width	420 mm

BALANCE F
ART. NO. S3.10.13.501.0

BALANCE F&B (with brake)
ART. NO. S3.11.13.501.0



BALANCE PLUS F

The BALANCE PLUS features adjustable armrests for optimal comfort and hygiene.

- Height and length of armrests can be customized to prevent tired arms
- Made of high-quality manmade material, the armrests are easy to clean and maintain high standards of hygiene
- Plain black design complements any color version of the chair

BALANCE PLUS F
ART. NO. S3.10.13.511.0

BALANCE PLUS F&B (with brake)
ART. NO. S3.11.13.511.0



BALANCE SUPREME F

The BALANCE SUPREME is equipped with premium armrests that are padded and multifunctional.

- The armrests are adjustable in height, length, and can be moved to the side, providing a comfortable and customizable sitting experience
- Select the color of the armrest cushions to match the seat surface

BALANCE SUPREME F
ART. NO. S3.10.13.521.0

BALANCE SUPREME F & B (with brake)
ART. NO. S3.11.13.521.0

Microsurgery chair **BALANCE ADVANCE**

The comfortable and ergonomic solution for precision work

Crafted in close collaboration with microsurgeons, the BALANCE ADVANCE operating chair is designed to provide maximum comfort and support during microsurgery procedures. With its ergonomic design and optimal working position, this chair enables you to work with precision and ease.



Available with and without brake

BALANCE ADVANCE F
ITEM No. S3.10.13.5310

BALANCE ADVANCE F&B (with brake)
ITEM No. S3.11.13.5310



Experience unmatched comfort and ergonomics

- Three-directional adjustable armrests with lockable feature support shoulder and arm muscles
- Comfortable work position even during extended procedures
- Ergonomically shaped seat and back with patented seat inclination and high-quality cushions ensure a healthy sitting posture
- Easy height adjustment with foot release (gas spring)
- Comfortable chassis with antistatic wheels
- High-frequency seamless cushions and closed wheels for hygienic cleaning

Technical data

Min. height	555 mm
Max. height	745 mm
Chair base diameter	610 mm
Seat cushion length	460 mm
Seat cushion width	420 mm

Microsurgery chair **THRONUS T3**

A first-class seating experience

You spend most of your life in the operation room, which is why we believe that high-quality equipment is not a luxury but a necessity for your work. Designed in collaboration with surgeons, every detail of the THRONUS T3 has been perfected to ensure outstanding, ergonomically perfected seating comfort.

With the THRONUS T3, you can enjoy a comfortable, stable seating position that helps you focus on your work with maximum precision. Investing in the THRONUS T3 microsurgery chair means investing in your professional success and well-being. Experience the difference for yourself and always sit in first class with the THRONUS T3.



Intraoperative armrest adjustment

Memory function enables programming of two seat heights

Electromechanical height adjustment

Standout features

- Optimized sitting and working positions with adjustable seat and back sections
- Intraoperatively adjustable, stable armrests for quick and easy fine-tuning
- Electronic height adjustment with two memory positions, activated via foot release
- Ultimate hygiene with plastic-wrapped cushion and covered wheels
- Button-activated electrical braking system - with LED indicator for stable positioning.
- The solid and heavy design increases stability
- Replacement battery system with an average of 4 days' runtime for uninterrupted operation
- Premium-quality manufacturing for long-lasting durability

Technical data

Min. height	555 mm
Max. height	830 mm
Seat cushion length	400 mm
Seat cushion width	360 mm

THRONUS T3 ART. NO. T3.00.00



Base plate for foot switch- THRONUS "set"

- For storing external foot controls (e.g. microscope)
- Flexibly extendable forward and spreadable to the side

ART. NO. Z1.219

Operating stool **COMFORT**

The allrounder

Designed to enable medical staff to work with speed and agility, this maneuverable stool offers maximum flexibility for seamless use. With its practical foot-operated height adjustment, users can easily achieve the desired height while keeping their hands free.

Standout features

- Seamless seat cushions for enhanced comfort and easy cleaning, ensuring reliable hygiene
- Large, enclosed, antistatic wheels that facilitate fast and agile working, while maintaining optimal hygiene
- Optional locking brake for a secure and stable working position (F/B)
- Optional 100 mm column elevation for increased versatility and customizable height adjustment (ART. NO. E.001558)



Available with
COMFORT SWING
seat cushion



Height adjustment
with foot release,
ideal for sterile work



Technical data COMFORT F

Min. height	540 mm
Max. height	735 mm
Chair cross diameter	610 mm
Seat diameter	400 mm
Flexible height adjustment	510 - 710 mm

Available COMFORT models

COMFORT F

Height adjustment
by foot release

ART. NO. H3.10.11.100.0

COMFORT F/B

Height adjustment by foot
release, with brake

ART. NO. H3.11.11.100.0



Optionally
with brake

COMFORT SWING F

With COMFORT SWING seat
cushion, height adjustment by
foot release

ART. NO. H3.10.12.100.0

COMFORT SWING F/B

With COMFORT SWING seat
cushion, height adjustment by
foot release, with brake

ART. NO. H3.11.12.100.0



Height adjustment
with hand release

COMFORT H

Height adjustment by
hand release

ART. NO. H3.10.31.100.0

COMFORT SWING H

With COMFORT SWING seat
cushion, height adjustment via
hand release

ART. NO. H3.10.32.100.0

OR stool **COMFORT ASSIST**

Your new, comfortable support in the operating room.

The COMFORT ASSIST surgical stool is the ideal addition to any operating theatre or treatment room. The 360° swivelling backrest of the stool provides ergonomic support in a variety of sitting positions and enables quick and agile working.

Standout features

- 360° swivelling backrest for ergonomic support

Technical data

Min. height	540 mm
Max. height	680 mm
Chair cross diameter	610 mm
Seat diameter	400 mm

COMFORT ASSIST
ART. NO. H3.10.41.103.0





Our BRUMABA operating theatre chairs and stools at a glance



	COMFORT H	COMFORT F	COMFORT F&B	COMFORT ASSIST F
Article number	H3.10.31.100.0	H3.10.11.100.0	H3.11.11.100.0	H3.10.41.103.0
Chair cross ø mm	610	610	610	610
Chair height mm	475-680	540-735	540-735	475-680
Height adjustment	Hand release	Foot release	Foot release	Foot release
Brake			•	
Seat cushion round (mm)	400	400	400	400

Max. load 120 kg

COMFORT SWING	COMFORT SWING H	COMFORT SWING F	COMFORT SWING F&B	
Article number	H3.10.32.100.0	H3.10.12.100.0	H3.11.12.100.0	
Chair height mm	570-750	500-705	570-750	



	BALANCE F	BALANCE F&B
Article number	S3.10.13.501.0	S3.11.13.501.0
Chair cross ø mm	610	610
Chair height mm	500-695	500-695
Height adjustment	Foot release	Foot release
Brake		•
Seat cushion BALANCE mm	460x420	460x420

Max. load 120 kg



	BALANCE ADVANCE F	BALANCE ADVANCE F&B
Article number	S3.10.13.531.0	S3.11.13.531.0
Chair cross ø mm	610	610
Chair height mm	555-745	555-745
Height adjustment	Foot release	Foot release
Brake		•
Seat cushion BALANCE mm	460x420	460x420
Max. load	120 kg	
SLIM LINE	BALANCE ADVANCE SLIM F	BALANCE ADVANCE SLIM F&B
Article number	S3.10.23.641.0	S3.12.23.641.0
Seat cushion in mm	400x360	400x360
Chair height mm	490-630	490-630

	BALANCE SUPREME F	BALANCE SUPREME F&B
Article number	S3.10.13.521.0	S3.11.13.521.0
Chair cross ø mm	610	610
Chair height mm	520-715	520-715
Height adjustment	Foot release	Foot release
Brake		•
Seat cushion BALANCE mm	460x420	460x420
Max. load	120 kg	










Unleash your creativity


35 vibrant cushion colors

Our cushion colors are available for both operating tables and stools/chairs, allowing for seamless integration of your BRUMABA products into your operating room environment.



Scan QR code
for further
information





Color swatches may differ from the original color tone.



BRUMABA GmbH

Bürgermeister-Graf-Ring 17
D-82538 Geretsried, Germany

Phone: +49 (0) 8171 / 2672 - 0

Fax: +49 (0) 8171 / 2672 - 10

Email: info@brumaba.de

Subject to change without notice.

© BRUMABA GmbH
KA.0013-2025-12-11-EN

www.brumaba.com

MADE IN GERMANY



SINCE 1980

BRUMABA
OPERATING TABLE SYSTEMS