




General Catalogue Accessories

Your extras for more
functionality and comfort

ACCESSORIES

for operating tables



BRUMABA is a second-generation, family-owned company founded in 1980. **BRUMABA** develops, manufactures, and globally distributes operating tables, operating chairs, and patient transporters and appropriate accessories for specialty departments in hospitals, clinics, and doctor's offices.



BRUMABA'S CORE COMPETENCE: MOBILE, SPECIALIZED OPERATING TABLES.

BRUMABA is a professional and competent address as well as adviser for all matters concerning ergonomics and positioning on operating tables. We offer our customers products for specialized departments such as eye and oral surgery, OMS, and plastic surgery, dermatology, as well as for general operations, both for in-patients and out-patients.

Operating tables for operating without tiring and optimum patient positioning.

This complete catalog of accessories contains appropriate accessories for expanding and customizing your entire operating table system.



OVERVIEW OF ACCESSORIES

POSITIONING | HEAD

Head plate, large, 400x300 mm	8
Head plate, small, 320x240 mm	8
Horseshoe-shaped headrest (included as standard)	8
Horseshoe-shaped headrest, small	9
Horseshoe-shaped headrest with long clamping piece, 110 mm	9
Silicone head restraining belt for fastening on the headrest	9
Hair transplant support	10
Neck cushion, viscoelastic, black	10
Head plate "rubins pillow" with ball-and-socket joint	10
Headrest PRIMUS / VARIUS, detachable	11
Radiolucent headrest, detachable (surcharge)	11
Microsurgery ring 2.0, with integrated mounting bracket	12
Microsurgery wrist rest GENIUS EYE / EYE PRO	12
Lower arm support for Ophthalmology	13
Tongue holder, extended, with flat plate rail	13
TriFlex Assist	13

POSITIONING | ARM

Armrest with fastening belt	14
Armrest OP - PRO, with fastening belt	14
Arm support lateral with fastening belt	14
Patient handle, hinged	15
Hand surgery table, with telescopic support	15
Telescopic support for hand surgery table	15
Back section cushion with arm support, visco	16
Side Guard, hinged	16
Horseshoe-shaped headrest for shoulder arthroscopy with head restraining belt	17
Back section extension for shoulder arthroscopy	17

POSITIONING | TORSO

Lateral positioning support, deep	18
Lateral positioning support, dorsal (back)	18
Body belt splitted, silicone with 2 mounting gliders	18
Leg holder PRO - OP (1 pair)	19

POSITIONING | LEG

Joint support	20
Arthroscopic leg support	20
Padded insert for arthroscopic surgery leg support	21
Standing plate (1 pair)	21

POSITIONING PADDINGS

Lordosis support padding	22
Half roll padding, small	22
Half roll padding, large	22
Stomach positioning padding, trapezoidal (1 piece)	23
Mattress, viscoelastic	23
Mattress, viscoelastic, short, with gynecological cutout	23

OPERATING TABLE | EXTENSIONS

Back section extension, narrow	24
Back section extension, wide	24
Back section extension, extra wide	25
Back section extension (GENIUS EYE/GENIUS EYE PRO)	25
Foot section extension Incl. Mounting bracket (GENIUS)	26
FOOT SECTION EXTENSION / SIDE EXTENSION NARROW (PRIMUS / VARIUS)	26
Side extension, radiolucent	26
Lower leg plate, combined (VARIUS)	27
Lower leg plate, combined, extended (VARIUS)	27
Lower leg plates, extended, 1 pair (VARIUS)	27

INSTRUMENT TABLES

Tray for motor, swifeling, on column	28
Instrument table, mobile	28
Dürr swivel slide on column	29
Instrument table on column	29
Mounting arm	29

ANESTHESIA / INFUSION STANDS

Anesthetic frame FLEX, foldable.....	30
Anesthetic frame, bendable – with oxygen supply.....	30
Anesthetic tube holder FLEX, foldable.....	31
Anesthetic frame, universal	32

INFUSION STANDS

Infusion stand with ball-and-socket joint	32
Infusion stand on the column, hinged.....	33
Infusion stand - swivel arm, extended (360 mm).....	33
Infusion stand mounted on instrument table	33

CONTROLS

Hand remote control type IV	34
Remote control holder on infusion stand	34
Remote control holder FLEX, foldable.....	35
Remote control holder with fixed clamp	35
Magnetic plate for remote control	35
Foot control, type IV.....	36
Foot control, type IV without memory.....	36
Platform for foot control type IV	37
Foot control, type IV fixed on undercarriage.....	37
Foot control, type IV EYE, without memory.....	37

MOUNTING | BRACKETS AND ADAPTERS

Mounting bracket.....	38
Mounting bracket special “XL”	38
Mounting glider with clamping	39
Mounting adapter, 25x10 mm.....	39
Mounting adapter, strengthened, 290 mm.....	39

MOUNTING | MOUNTING RAILS

Mounting rail, 250 mm head section.....	40
Mounting rail, 340 mm foot section (GENIUS).....	40
Mounting rail, 250mm massiv seat section (GENIUS).....	40
Mounting rail, 340 mm, solid back and foot section (GENIUS).....	41
Mounting rail, 125 mm Head section GENIUS EYE / EYE PRO).....	41
Mounting rails set (VENUS).....	41

TRANSPORT HANDLES

Transport handles, foldaway, head end, GENIUS (1 pair)	42
Transport handles, foldaway, head end PRIMUS, VARIUS (1 pair)	42
Transport handles, hinged on foot side (1 pair)	42
EasyCover (150 pieces – disposable, unsterile).....	43
Paper towel holder	43

HYGIENE / PRESERVATION / CARE

Protective cover for foot section, VISCO	44
Protective cover for foot section, VISCO, (VARIUS) (1 pair)	44

MAGNIFYING LAMPS

Magnifying lamp with swivel arm	45
---------------------------------------	----

ACCESSORIES ORDER

Wall mounting rail	46
Wall mount for hand surgery table.....	46
ACCESSORIES column.....	46
Accessories trolley.....	47
Transport trolley for Batteries	47

HEAD PLATE, LARGE, 400 x 300 mm

ART.-NO.: Z1.003



PADDING DIMENSIONS (LxWxH):
400x300x50 mm

CLAMPING PIECE HEIGHT: 110 mm


- For head positioning, lateral positioning
- Abdominal positioning in conscious patients
- Ball joint for seamless inclination, height-adjustable

Suitable for:

- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

HEAD PLATE, SMALL, 320 x 240 mm

ART.-NO.: Z1.109



PADDING DIMENSIONS (LxWxH):
320x240x50 mm

CLAMPING PIECE HEIGHT: 110 mm

- For head positioning, lateral positioning
- Abdominal positioning in conscious patients
- Ball joint for seamless inclination, height-adjustable

Suitable for:

- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

HORSESHOE-SHAPED HEADREST (INCLUDED AS STANDARD)

ART.-NO.: Z1.043



PADDING DIMENSIONS (LxWxH):
250x240x110 mm

CLAMPING PIECE HEIGHT: 65 mm


- For a defined and stable head position
- Ball-and-socket joint, height adjustable

Suitable for:

- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

HORSESHOE-SHAPED HEADREST, SMALL

ART.-NO.: Z1.043A



PADDING DIMENSIONS (LxWxH):
200x210x80 mm

CLAMPING PIECE HEIGHT: 65 mm

- For good access to the ears
- For head positioning in children
- Ball-and-socket hinge, height adjustable,

Suitable for:

- » PRIMUS » GENIUS
- » VARIUS » GENIUS EYE
- » GENIUS EYE PRO

HORSESHOE-SHAPED HEADREST WITH LONG CLAMPING PIECE, 110 mm

ART.-NO.: Z1.249



CLAMPING PIECE HEIGHT: 110 mm

- For Bekhterev's disease patients
- Defined and stable head position
- Ball-and-socket hinge, height adjustable

Suitable for:

- » PRIMUS » GENIUS
- » VARIUS » GENIUS EYE
- » GENIUS EYE PRO

SILICONE HEAD RESTRAINING BELT FOR FASTENING ON THE HEADREST

ART.-NO.: Z1.189



BELT SIZE (LxW):
900x50 mm

- Hygienic silicone belt for the fixation of the patient's head

Suitable for:

- » PRIMUS » GENIUS
- » VARIUS » GENIUS EYE
- » GENIUS EYE PRO

HAIR TRANSPLANT SUPPORT

ART.-NO.: Z1.032



PADDING DIMENSIONS (BxL):
90x230 mm

- For optimum access to the upper and rear head area
- Ball-and-socket hinge, height adjustable

Suitable for:

- » PRIMUS » GENIUS
- » VARIUS » GENIUS EYE
- » GENIUS EYE PRO

NECK CUSHION, VISCOELASTIC, BLACK

ART.-NO.: Z1.250



PADDING DIMENSIONS (BxL):
200 x 200 mm

CLAMPING PIECE HEIGHT: 65 mm

- For extra soft neck positioning
- Ball-and-socket hinge, height adjustable

Suitable for:

- » PRIMUS » GENIUS
- » VARIUS » GENIUS EYE
- » GENIUS EYE PRO

HEAD PLATE "RUBINS PILLOW" WITH BALL-AND-SOCKET JOINT

ART.-NO.: X.000184



PADDING DIMENSIONS (LxWxH):
470 x 250 x 100 mm

CLAMPING PIECE HEIGHT: 110 mm

- For defined head positioning
- Armrest for the surgeon during eye surgery in 12 o'clock position
- Ball-and-socket hinge, height adjustable,

Suitable for:

- » PRIMUS » GENIUS
- » VARIUS » GENIUS EYE
- » GENIUS EYE PRO

HEADREST PRIMUS / VARIUS, DETACHABLE

ART.-NO.: X.000209



- For a quick change of the head section
- Good access for the anesthetist, e.g. during knee or gynecological surgery
- Solid stainless steel construction

Suitable for:

- » VARIUS
- » PRIMUS

RADIOLUCENT HEADREST, DETACHABLE (SURCHARGE) **ART.-NO.: X.000235**




OPTIONAL: PADDING FOR
HEAD SECTION RADIO LUCENT
(AS SHOWN)

ORDER-NO.: E.002553

- For a quick change of the head section
- Good for intra-operative X-ray, e.g. ENT, OMS

Suitable for:

- » VARIUS
- » GENIUS
- » PRIMUS



**MICROSURGERY RING 2.0,
WITH INTEGRATED MOUNTING BRACKET**

ART.-NO.: Z1.240



DIMENSIONS:

INTERIOR B: 325 mm

EXTERIOR B: 440 mm

PADDING WIDTH: 40 mm


- For precise operating and operating without tiring
- Ideal for microscope surgery
- Ball-and-socket hinge, height adjustable
- Soft padding

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

MICROSURGERY WRIST REST GENIUS EYE / EYE PRO

ART.-NO.: Z1.257



DIMENSIONS:

INTERIOR B: 292 mm

EXTERIOR B: 372 mm


- For precise operating and operating without tiring
- Ideal for microscope surgery
- Ball-and-socket hinge, height adjustable

Suitable for:

- » GENIUS EYE
- » GENIUS EYE PRO

LOWER ARM SUPPORT FOR OPHTHALMOLOGY

ART.-NO.: Z1.298 | GENIUS EYE PRO



- Ergonomic lower armrest for the surgeon during eye surgery in 12 o'clock position
- For operating without tiring during ophthalmological procedures
- Soft padding

Suitable for:

» GENIUS EYE PRO

TONGUE HOLDER, EXTENDED, WITH FLAT PLATE RAIL

ART.-NO.: Z1.061



DIMENSIONS (W x H):
460 x 390 – 670 mm

OPTIONAL: TONGUE HOLDER, SHORT, WITH RAISED PLATE RIM (AS SHOWN)

ORDER-NO.: Z1.061

- For ENT
- Large individual range by way of two ball-and-socket hinges
- Chain for tongue and retractor holder
- instrument tray
- Height-adjustable

Suitable for:

» VARIUS » GENIUS
» PRIMUS

TRIFLEX ASSIST

ART.-NO.: Z1.363



DIMENSIONS:
47cm length, unfolded
max. ø 4.5cm

- Ideal for single practitioners
- Magnetic holder
- Multidirectional movement 360° ball joint
- Only available in combination with microsurgical ring

Suitable for:

» VARIUS » GENIUS
» PRIMUS

ARMREST WITH FASTENING BELT



ART.-NO.: Z1.001

PADDING DIMENSIONS (LxWxH):
490x140x40 mm

RESTRAINING BELT:
650x50 mm

- For ergonomic arm positioning
- Eccentric lever operates two ball-and-socket hinges, for fast, seamless adjustment
- Also for prone position
- Integrated impact protection
- Detachable padding, Silicone restraining belt

Suitable for:

» all BRUMABA Operating Tables

ARMREST OP - PRO, WITH FASTENING BELT



PADDING DIMENSIONS (LxWxH):
490X140X50 mm

RESTRAINING BELT:
650X50 mm

- Simple length adjustmentFor ergonomic arm positioning
- Eccentric lever operates two ball joints for quick, stepless adjustment
- Also for prone position
- Integrated impact protection
- Removable padding, silicone fixation strap, suitable for all BRUMABA operating tables

Suitable for:

» all BRUMABA Operating Tables

ARM SUPPORT LATERAL WITH FASTENING BELT



ART.-NO.: Z1.002

PADDING DIMENSIONS (LxWxH):
240x140x40 mm

- For slender arm apposition
- Folds down by way of swivel joint
- Detachable padding, Silicone restraining belt

Suitable for:

» all BRUMABA Operating Tables

PATIENT HANDLE, HINGED



PATIENT HANDLE FOR "VARIUS, PRIMUS, GENIUS"

ART.-NO.: Z1.281

PATIENT HANDLE FOR "VENUS"

ART.-NO.: Z1.297

■ For keeping the patient in the prone position

Suitable for:

» VARIUS » GENIUS
» PRIMUS » VENUS

HAND SURGERY TABLE, WITH TELESCOPIC SUPPORT



PADDING DIMENSIONS (LxWxH):
710 X 420 mm

■ For ergonomic positioning of the patient's arm
■ Ball joint
■ X-ray compatible
■ Slim design for maximum legroom

Suitable for:

» VARIUS » VENUS
» PRIMUS » GENIUS

TELESCOPIC SUPPORT FOR HAND SURGERY TABLE



HEIGHT ADJUSTABLE FROM::
570-1030 mm

■ For higher stability requirements on the hand table
■ Height adjustable

Suitable for:

» VARIUS » VENUS
» PRIMUS » GENIUS

BACK SECTION CUSHION WITH ARM SUPPORT, VISCO

ART.-NO.: Z1.259



PADDING WIDTH:
720 mm

■ For simple arm positioning

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

SIDE GUARD, HINGED

ART.-NO.: Z1.187



DIMENSIONS (W x H):
350 x 310 mm

■ For additional safety during transport of the patient

■ Folds down by way of swivel joint

Suitable for:

- » all BRUMABA Operating Tables

HORSESHOE-SHAPED HEADREST FOR SHOULDER ARTHROSCOPY WITH HEAD RESTRAINING BELT

ART.-NO.: Z1.154



PADDING DIMENSIONS (LxWxH):
200x210x80 mm

CLAMPING PIECE HEIGHT: 110 mm

- For laterally offset positioning of the head
- For free access to the left and right shoulder area
- In combination with the head cradle, the horseshoe-shaped headrest adapts automatically during adjustment
- Ball-and-socket hinge, height adjustable, silicone restraining belt

Suitable for:

- » VARIUS
- » PRIMUS

BACK SECTION EXTENSION FOR SHOULDER ARTHROSCOPY

ART.-NO.: Z1.021



PADDING DIMENSIONS (LxW x H):
270x150x60 mm

- For shoulder operations left and right
- Easy assembly

Suitable for:

- » VARIUS
- » PRIMUS

LATERAL POSITIONING SUPPORT, DEEP

ART.-NO.: Z1.110



PADDING DIMENSIONS (LxWxH):
180x300 x 50 mm

CLAMPING PIECE HEIGHT: **120 mm**

OPTIONAL:
LATERAL POSITIONING SUPPORT
ORDER-NO.: Z1.054

CLAMPING PIECE HEIGHT: **65 mm**

- For stable lateral positioning of the patient
- Ball-and-socket hinge for seamless adjustment

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

LATERAL POSITIONING SUPPORT, DORSAL (BACK)

ART.-NO.: Z1.152



PADDING DIMENSIONS (LxWxH):
240x140 x 40 mm

HEIGHT-ADJUSTABILITY **120 mm**

- For slender support of the patient in lateral position
- Two joints for individual adaptation
- Height-adjustable

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

BODY BELT SPLITTED, SILICONE WITH 2 MOUNTING GLIDERS

ART.-NO.: Z1.112



DIMENSION (LxW):
1440x80 mm

- For securing the patient
- Hygienic silicon material
- Reduced pressure load

Suitable for:

- » all BRUMABA Operating Tables

LEG HOLDER PRO - OP (1 PAIR)

ART.-NO.: Z1.366



PADDING DIMENSIONS (LxWxH):

220x330X90 mm

RESTRAINING BELT DIMENSIONS

(LxB): **650x50 mm**


HEIGHT-ADJUSTABILITY: **130 mm**

- For ergonomic leg positioning
- Three-dimensional mobility
- Easy to handle by just one person
- Rocker function for automatic adaptation of the leg band against positional damage
- Ball-and-socket joints, height adjustable, silicon restraining belt

Suitable for:

- » VARIUS » VENUS
- » PRIMUS » GENIUS

POSITIONING | SURGICAL LEG / KNEE



JOINT SUPPORT

ART.-NO.: Z1.028



PADDING DIMENSIONS (BxL):
90x230 mm

HEIGHT ADJUSTABLE FROM:
380–480 mm

- For surgery on the foot or shin
- Ball-and-socket joint
- Height-adjustable

Suitable for:

- » VARIUS » VENUS
- » PRIMUS » GENIUS

ARTHROSCOPIC LEG SUPPORT

ART.-NO.: Z1.016



MAX. B THIGH: 270 mm

INCL. MOUNTING BRACKET XL

ORDER-NO.: Z1.088


INCL. MOUNTING ADAPTER

ORDER-NO.: Z1.047

- for knee arthroscopic surgery
- Suitable for left and right knee
- Simple thigh adaptation
- Also for cruciate ligament surgery
- Can be rotated 360°
- Height-adjustable

Suitable for:

- » VARIUS



**PADDED INSERT FOR ARTHROSCOPIC SURGERY
LEG SUPPORT**

ART.-NO.: Z1.155

PADDING THICKNESS: 25 mm

■ For knee arthroscopic surgery
without tourniquet

Suitable for:
» VARIUS



STANDING PLATE (1 PAIR)

ART.-NO.: Z1.057

PADDING DIMENSIONS (B x L):
180 x 300 mm

■ For procedures in the reverse
Trendelenburg position, e.g. phlebology

Suitable for:
» VARIUS » VENUS

FOR „PRIMUS“

ART.-NO.: X.000359



POSITIONING PADDINGS

LORDOSIS SUPPORT PADDING

ART.-NO.: Z1.140



PADDING DIMENSIONS (LxWxH):
400x210x60 mm

- Soft material and hygienic, seamless cover
- For the support of the lordosis

Suitable for:
» all BRUMABA Operating Tables

HALF ROLL PADDING, SMALL

ART.-NO.: Z1.100



PADDING DIMENSIONS (LxWxH):
400x140x70 mm

- Soft material and hygienic, seamless cover
- Comfortable shape for various purposes

Suitable for:
» all BRUMABA Operating Tables

HALF ROLL PADDING, LARGE

ART.-NO.: Z1.101



PADDING DIMENSIONS (LxWxH):
555x210x120 mm

- Soft material and hygienic, seamless cover
- Various purposes

Suitable for:
» all BRUMABA Operating Tables

POSITIONING PADDINGS

STOMACH POSITIONING PADDING, TRAPEZOIDAL (1 PIECE)

ART.-NO.: Z1.098



PADDING DIMENSION (LxW):
220x300 mm

MAX. DEPTH: 110 mm

- Soft material and hygienic, seamless cover
- For ergonomic stomach positioning to relieve the abdominal area

Suitable for:

» all BRUMABA Operating Tables

MATTRESS, VISCOELASTIC

ART.-NO.: Z1.171



PADDING DIMENSIONS (LxWxH):
1660x570x40 mm

- For long operations to prevent bedsores and for extra comfort
- Anti-static

Suitable for:

» all BRUMABA Operating Tables

MATTRESS, VISCOELASTIC, SHORT, WITH GYNECOLOGICAL CUTOUT

ART.-NO.: Z1.213




PADDING DIMENSIONS (LxWxH):
950 x 550 mm

- For long operations to prevent bedsores and for extra comfort
- Anti-static

Suitable for:

» all BRUMABA Operating Tables



BACK SECTION EXTENSION, NARROW

**ART.-NO.: Z1.006 | PRIMUS / VARIUS
Z1.183 | GENIUS**



PADDING DIMENSION (H x W):
100 x 420 mm

- For correct positioning of the shoulders in large patients
- Included as standard with PRIMUS, VARIUS and GENIUS

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

BACK SECTION EXTENSION, WIDE

**ART.-NO.: Z1.005 | PRIMUS / VARIUS
Z1.184 | GENIUS**



PADDING DIMENSION (H x W):
150 x 420 mm

- For correct positioning of the shoulders in large patients

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

BACK SECTION EXTENSION, EXTRA WIDE

**ART.-NO.: Z1.212 | PRIMUS / VARIUS
Z1.272 | GENIUS**



PADDING DIMENSION (HxW):
150 x 550 mm

- For correct positioning of the shoulders in large patients
- Extra wide

Suitable for:

- » VARIUS » GENIUS
- » PRIMUS

BACK SECTION EXTENSION (GENIUS EYE/GENIUS EYE PRO)

ART.-NO.: Z1.279



PADDING DIMENSION (HxW):
150 x 530 mm

- For correct positioning of the shoulders in large patients

Suitable for:

- » GENIUS EYE
- » GENIUS EYE PRO

OPERATING TABLE | EXTENSIONS

FOOT SECTION EXTENSION INCL. MOUNTING BRACKET (GENIUS)

ART.-NO.: Z1.291 | NARROW
ART.-NO.: Z1.286 | WIDE



PADDING DIMENSION (HxW):
 NARROW 100 x 420 mm
 WIDE 150 x 420 mm

- For particularly large patients
- With edge protection plate
- Additional mounting bracket required

Suitable for:

» GENIUS

Side extension for GENIUS:

ART.-NR.: Z1.263 | NARROW

ART.-NR.: Z1.264 | WIDE

FOOT SECTION EXTENSION / SIDE EXTENSION NARROW (PRIMUS / VARIUS)

ART.-NO.: Z1.321 | NARROW*
ART.-NO.: Z1.268 | WIDE



PADDING DIMENSION (HxW): :
 NARROW 100 x 420 mm
 WIDE 150 x 420 mm

- For particularly large patients or widening the lying surface
- With edge protection plateAdditional mounting bracket required

Suitable for:

» VARIUS » PRIMUS

Without mounting adapter

*with mounting adapter

SIDE EXTENSION, RADIOLUCENT

ART.-NO.: Z1.055



DIMENSION (LxW x H):
 410 x 225 X 55 mm

- For evacuation of the extremities during operations and X-rays
- Additional mounting bracket required

Suitable for:

» VARIUS » VENUS
 » PRIMUS » GENIUS

LOWER LEG PLATE, COMBINED (VARIUS)

ART.-NO.: Z1.170



PADDING DIMENSIONS | VISKO
(LxWxH):
500x555x60 mm

- For use during inverse positioning
e.g. with external mounting
ORDER-NO.: Z1.194 for a
homogeneous X-ray image
- Connecting plate for the two-part foot
section

Suitable for:

» VARIUS

LOWER LEG PLATE, COMBINED, EXTENDED (VARIUS)

ART.-NO.: Z1.169



PADDING DIMENSIONS | VISKO
(LxWxH):
800x555x60 mm

- For use during inverse positioning
- Extension of the X-ray range
- For a homogeneous X-ray image

Suitable for:

» VARIUS

LOWER LEG PLATES, EXTENDED, 1 PAIR (VARIUS)

ART.-NO.: Z1.020




PADDING DIMENSIONS (LxWxH):
800x270x50 mm

- Extension of the leg plates e.g. after
lithotomy position or knee surgery

Suitable for:

» VARIUS

MOTOR MOUNTING BRACKET / INSTRUMENT TABLES



TRAY FOR MOTOR, SWIFELING, ON COLUMN

ART.-NO.: X.000040



PLATE SIZE: 400x360 mm

■ For storing the surgical or piezosurgical units

Suitable for:

» all BRUMABA Operating Tables

INSTRUMENT TABLE, MOBILE

ART.-NO.: Z1.172



PLATE SIZE: 600x400 mm

OPTIONAL:

ART.-NO.: Z1.282

PLATE SIZE: 700x500 mm

HEIGHT ADJUSTABLE FROM:
860–1100 mm

- For the storage of instruments within easy reach
- Mobile with hygienic, smooth-running double castors
- One-handed height adjustment via gas spring
- Brake against rolling away

INSTRUMENT TABLES

DÜRR SWIVEL SLIDE ON COLUMN

ART. NO. Z1.372



DÜRR SWIVEL SLIDE WITH INSTRUMENT TABLE ON COLUMN

PLATE SIZE: 600x400 mm

ART.-NO.: Z1.373

PLATE SIZE: 700x500 mm

ART.-NO.: Z1.374

- Swiveling and extendable by approx. 310 mm
- Includes additional swivel joint
- Total length 700 mm
- WITHOUT hose reel

Suitable for:

» PRIMUS » GENIUS

INSTRUMENT TABLE ON COLUMN

ART.-NO.: Z1.275



PLATE SIZE: 600x400 mm

OPTIONAL:

ART.-NO.: Z1.278

PLATE SIZE: 700x500 mm

HEIGHT-ADJUSTABLE ABOVE TABLE LEVEL: 300–540 mm

- For the storage of instruments within ergonomic reach above the patient
- One-handed adjustment via gas spring
- Plate 360° rotatable, arm 180° tiltable
- Left and right mounting possible

Suitable for:

» all BRUMABA Operating Tables

MOUNTING ARM



DOUBLE SUPPORT ARM ON BOTH SIDES

FOR INSTRUMENT TABLE AND INFUSION STAND

ART.-NO.: M.000380


DOUBLE SUPPORT ARM OPERATING TABLE INSTRUMENT TABLE COMPLETE:

ART.-NO.: M.000843

Suitable for:

» VARIUS » GENIUS
» PRIMUS

ANESTHESIA / INFUSION STANDS



ANESTHETIC FRAME FLEX, FOLDABLE

ART.-NO.: Z1.273



LENGTH: 960 mm

- For the supply of oxygen to the patient's head
- External hygienic tube guide on clips for 6 mm and 18 mm tubes (recommendation)
- Spacer for surgical drape
- Flexible gooseneck for individual adjustment
- Folds down by way of swivel joint

Suitable for:

» all BRUMABA Operating Tables

ANESTHETIC FRAME, BENDABLE – WITH OXYGEN SUPPLY

ART.-NO.: Z1.248



LENGTH: 1010 mm

- For the supply of oxygen to the patient's head
- Oxygen tube connection at gooseneck with outlet at head
- Spacer for surgical drape
- Flexible gooseneck for individual adjustment



Suitable for:

» all BRUMABA Operating Tables

ANESTHESIA / INFUSION STANDS

ANESTHETIC TUBE HOLDER FLEX, FOLDABLE

ART.-NO.: Z1.334




I For anesthetic tube holder

Suitable for:

» all BRUMABA Operating Tables

ANESTHESIA / INFUSION STANDS



ANESTHETIC FRAME, UNIVERSAL

ART.-NO.: Z1.015



WIDTH: 680 – 1130 mm
HEIGHT-ADJUSTABLE: 270 mm

- For all standard covering systems in surgery
- Length and height adjustable
- Ball-and-socket hinge
- Hook for infusion bottle

Suitable for:

» all BRUMABA Operating Tables

INFUSION STAND WITH BALL-AND-SOCKET JOINT

ART.-NO.: Z1.038



HEIGHT: 1110 – 1680 mm
WITH 4 HOOKS
MAX. PERMISSIBLE LOAD:
2 KG PER HOOK

- Mounting on rail system
- One-handed height adjustment
- Ball-and-socket hinge
- 4 bottle holders

Suitable for:

» all BRUMABA Operating Tables

INFUSION STANDS

INFUSION STAND ON THE COLUMN, HINGED

ART.-NO.: Z1.277



HEIGHT: 1110–1680 mm

WITH 4 HOOKS

MAX. PERMISSIBLE LOAD:
2 KG PER HOOK

- Mounting on the column
- Independent of the angle of inclination setting of the patient surface
- One-handed height adjustment
- Tilttable
- 4 bottle holders

Suitable for:

» all BRUMABA Operating Tables

INFUSION STAND - SWIVEL ARM, EXTENDED (360 mm)

ART.-NO.: X.000214



OPTIONAL TO **ORDER-NO.: Z1.181**

LENGTH: 360 mm

- For a larger tilting range of the infusion stand
- ORDER-NO.: Z1.181**

Suitable for:

» all BRUMABA Operating Tables

INFUSION STAND MOUNTED ON INSTRUMENT TABLE

ART.-NO.: Z1.280




- Hinged on instrument table on column
- ORDER-NO.: Z1.275 / Z1.278**
- One-handed height adjustment
 - 2 bottle holders

Suitable for:

» all BRUMABA Operating Tables

CONTROLS



HAND REMOTE CONTROL TYPE IV

ART.-NO.: Z1.118




- Included as standard with delivery of the operating table
- Magnetic back
- Carry out all individual settings
- Retrieval of all 8 memory positions

Suitable for:

- » PRIMUS
- » GENIUS
- » VARIUS
- » GENIUS EYE
- » GENIUS EYE PRO

REMOTE CONTROL HOLDER ON INFUSION STAND

ART.-NO.: Z1.243



- For attaching the hand-held remote control to the infusion stand

Suitable for:

- » all BRUMABA Operating Tables

REMOTE CONTROL HOLDER FLEX, FOLDABLE

ART.-NO.: Z1.313



GOOSENECK LENGTH : 300 mm


■ Remote control holder FLEX, foldable

Suitable for:

» all BRUMABA Operating Tables

REMOTE CONTROL HOLDER WITH FIXED CLAMP

ART.-NO.: Z1.307



■ Mounting on rail system underneath patient surface

Suitable for:

» all BRUMABA Operating Tables

MAGNETIC PLATE FOR REMOTE CONTROL

ART.-NO.: Z1.045




■ Included as standard with delivery of every operating table

■ Attachment by way of rear adhesive surface

Suitable for:

» all BRUMABA Operating Tables

CONTROLS



FOOT CONTROL, TYPE IV

ART.-NO.: Z1.185
ART.-NO.: Z1.303 | VENUS

WEIGHT: 4.25 kg



- Carry out all individual settings
- Retrieval of all 8 memory positions
- Lifting bar for lifting or moving the foot control with the foot
- Can be operated in sterile condition

Suitable for:

- » PRIMUS
- » GENIUS
- » VARIUS
- » GENIUS EYE
- » GENIUS EYE PRO

FOOT CONTROL, TYPE IV WITHOUT MEMORY

ART.-NO.: Z1.271

WEIGHT: 4.25 kg



- Carry out all individual settings
- Lifting bar for lifting or moving the foot control with the foot
- No retrievable memory positions

Suitable for:

- » all BRUMABA Operating Tables

PLATFORM FOR FOOT CONTROL TYPE IV

ART.-NO.: Z1.274



- For positioning the foot control on the undercarriage
- 4-EASY-DRIVE, 5-WHEEL-DRIVE

Suitable for:

- » all BRUMABA Operating Tables

FOOT CONTROL, TYPE IV FIXED ON UNDERCARRIAGE

ART.-NO.: Z1.225



WEIGHT: 4.25 kg

- Carry out all individual settings
- Retrieval of all 8 memory positions
- Foot control fixed on undercarriage
- 4-EASY-DRIVE, 5-WHEEL-DRIVE

Suitable for:

- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

FOOT CONTROL, TYPE IV EYE, WITHOUT MEMORY

ART.-NO.: Z1.226




WEIGHT: 4.25 kg

- Adjustment of head section, back section and height possible
- Foot control fixed on undercarriage
- No retrievable memory positions

Suitable for:

- » all BRUMABA Operating Tables



MOUNTING BRACKET

ART.-NO.: Z1.004



FLANGE DIMENSION (LxB):
30x10 mm

- For fastening of BRUMABA accessories
- Safe 4-point clamping

Suitable for:
» all BRUMABA Operating Tables

MOUNTING BRACKET SPECIAL "XL"

ART.-NO.: Z1.088



FLANGE DIMENSION (LxB):
30x15 mm

- For fastening of BRUMABA accessories
- Safe 4-point clamping

Suitable for:
» all BRUMABA Operating Tables

MOUNTING GLIDER WITH CLAMPING

ART.-NO.: Z1.093



FLANGE DIMENSION (LxB):
30x10 mm


■ For the permanent fastening of accessories

Suitable for:

» all BRUMABA Operating Tables

MOUNTING ADAPTER, 25x10 mm

ART.-NO.: Z1.047



DIMENSIONS (LxHxB):
200x25x10 mm

OPTIONAL: MOUNTING ADAPTER,
30x10 mm (SPECIAL)

ART.-NO.: Z1.047A

DIMENSIONS (LxHxB):
200x30x10 mm


■ For fastening accessories of other manufacturers (standard rails)

Suitable for:

» all BRUMABA Operating Tables

MOUNTING ADAPTER, STRENGTHENED, 290 mm

ART.-NO.: Z1.341



DIMENSIONS (LxHxB):
290x25x10 mm

■ For fastening accessories of other manufacturers (standard rails)

Suitable for:

» all BRUMABA Operating Tables

MOUNTING RAIL, 250 mm HEAD SECTION

ART.-NO.: Z1.151



■ For head section left

Suitable for:

- » GENIUS
- » VENUS

MOUNTING RAIL, 340 mm FOOT SECTION (GENIUS)

ART.-NO.: Z1.079




■ For foot section left or right

Suitable for:

- » GENIUS

MOUNTING RAIL, 250MM MASSIV SEAT SECTION (GENIUS)

ART.-NO.: X.000087



■ For seat section, to attach leg holders
ART.-NO.: Z1.087

Suitable for:

- » GENIUS

MOUNTING RAIL, 340 MM, SOLID BACK AND FOOT SECTION (GENIUS)

ART.-NO.: X.000036




■ For back and foot section

Suitable for:

» GENIUS

MOUNTING RAIL, 125 mm HEAD SECTION (GENIUS EYE / EYE PRO)

ART.-NO.: Z1.262



■ E.g. for mounting the
microsurgery ring GENIUS EYE
ART.-NO.: Z1.257


Suitable for:

» GENIUS EYE

» GENIUS EYE PRO

MOUNTING RAILS SET (VENUS)

ART.-NO.: Z1.304



■ For back, seat and foot section

Suitable for:

» VENUS

TRANSPORT HANDLES

TRANSPORT HANDLES, FOLDAWAY, HEAD END, GENIUS (1 PAIR)

ART.-NO.: Z1.210 GENIUS




! For the back-friendly moving of the operating table in chair position

Suitable for:

- » GENIUS
- » GENIUS EYE PRO

TRANSPORT HANDLES, FOLDAWAY, HEAD END PRIMUS, VARIUS (1 PAIR)

ART.-NO.: Z1.211



! For the back-friendly moving of the operating table in chair position

Suitable for:

- » VARIUS
- » PRIMUS

TRANSPORT HANDLES, HINGED ON FOOT SIDE (1 PAIR)

ART.-NO.: Z1.164




! Für das rückenschonende Schieben des OP-Tisches von der Fußseite

Geeignet für alle BRUMABA OP-Tische

MONTAGEKLOBEN XL
ERFORDERLICH

ART.-NO.: Z1.088

HYGIENE / PRESERVATION / CARE



EASYCOVER (150 PIECES – DISPOSABLE, UNSTERILE)

ART.-NO.: Z1.138



DIMENSION (LxW):
2000x720 mm

PACKAGING UNIT: 150 PIECES
COLOR: WHITE

- Watertight disposable surgery bed drape for maximum hygiene
- Head pocket prevents cover from slipping
- Paper coating for comfort while lying

Suitable for:

» all BRUMABA Operating Tables

PAPER TOWEL HOLDER

VARIUS Art.-No.: Z1.294

VENUS Art.-No.: Z1.324

PRIMUS/GENIUS Art.-No.: Z1.051



DIMENSIONS (W):
560 mm

USE OF PAPER TOWELS UP TO
MAX. WIDTH OF 560 mm

- For attaching paper towels under the foot section

Suitable for:

- » PRIMUS
- » VARIUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

PROTECTIVE COVER FOR FOOT SECTION, VISCO

ART.-NO.: Z1.332



■ For the mechanical protection of the foot section padding e.g. outdoor footwear


■ Easy disassembly for cleaning

Suitable for:

- » PRIMUS
- » VENUS
- » GENIUS
- » GENIUS EYE
- » GENIUS EYE PRO

PROTECTIVE COVER FOR FOOT SECTION, VISCO, (VARIUS) (1 PAIR)

ART.-NO.: Z1.148A



■ For the mechanical protection of the foot section padding e.g. outdoor footwear

■ Easy disassembly for cleaning

Suitable for:

- » VARIUS

MAGNIFYING LAMPS

MAGNIFYING LAMP WITH SWIVEL ARM

ART.-NO.: Z1.217



LENS DIMENSION B: 160 mm

LENS DIOPTR 3.5

COLOR RENDITION: 6500 K

LED ENERGY EFFICIENCY: A+

Suitable for:

- » PRIMUS
- » GENIUS
- » VARIUS
- » GENIUS EYE PRO

Optional:

ART.-NO.: Z1.300 | VENUS

ACCESSORIES ORDER

WALL MOUNTING RAIL

ART.-NO.: Z1.162



LENGTH: 100 mm

- For fastening accessories to the wall
- 5 gliders included

Suitable for:

» all BRUMABA Operating Table accessories

GLIDERS FOR WALL MOUNTING /
ACCESSORIES TROLLEY:

ART.-NO.: Z1.161

ACCESSORIES COLUMN

ART.-NO.: Z1.117



HEIGHT: 1080 mm

B BASE: 620 mm

- For storage of up to 8 accessories
- Mobile

Suitable for:

» all BRUMABA Operating Table accessories

Optional accessory:

STORAGE PLATE **ART.-NO.: Z1.153**

ACCESSORIES ORDER

ACCESSORIES TROLLEY

ART.-NO.: Z1.060



DIMENSIONS (LxW x H):
730x550x1000 mm

RAIL LENGTH: 725 mm

- For the storage of accessories
- Mobile
- 12 gliders included
- Extension for hand table and arthroscopic surgery leg support

Suitable for:

» all BRUMABA Operating Table accessories

GLIDERS FOR WALL MOUNTING /
ACCESSORIES TROLLEY:

ART.-NO.: Z1.161

TRANSPORT TROLLEY FOR BATTERIES

ART.-NO.: Z1.175



DIMENSIONS (LxW x H):
300x450x800 mm

- For easy transportation of two battery packs
 - Holder for the charger
- ART.-NO.: Z1.134**

Suitable for:

» VARIUS
» PRIMUS

BRUMABA GmbH

Bürgermeister-Graf-Ring 17
D-82538 Geretsried, Germany

Phone: +49 (0) 8171 / 2672 - 0

Fax: +49 (0) 8171 / 2672 - 10

Email: info@brumaba.de

Subject to change without notice.

© BRUMABA GmbH

KA.007-2026-06-23-EN

www.brumaba.com

MADE IN GERMANY



SINCE 1980


BRUMABA
ERGONOMIC MEDICAL SOLUTIONS